

TO: WMI Authorized Dealers

FROM: White Mountain Imaging, Technical Support

DATE: 02.12.2018

SUBJECT: Permanently Deleting Images from Image Suite


Product: Image Suite

Resolution:

1. Open IS Web PACS
2. Select the patient(s) to be deleted by putting a check in the box next to the study

Select All			Patient ID	Name
<input checked="" type="checkbox"/>				
<input checked="" type="checkbox"/>	+		83472	Brown,Julie
<input checked="" type="checkbox"/>	+		43344	booth,ronald ^{AA}
<input checked="" type="checkbox"/>	+		83054	GLENN,STEVEN ^{AA}
<input checked="" type="checkbox"/>	+		74297	Clark,Nancy ^{AA}
<input checked="" type="checkbox"/>	+		23967	HONCOOP,WILLIAM ^{AA}
<input checked="" type="checkbox"/>	+		83490	Hungerfod,Richard ^{AA}
<input checked="" type="checkbox"/>	+		37274	homann,peter ^{AA}
<input checked="" type="checkbox"/>	+		25908	davies,monica ^{AA}
<input checked="" type="checkbox"/>	+		80065	BOWEN,BRUCE ^{AA}
<input checked="" type="checkbox"/>	+		37274	homann,peter ^{AA}
<input checked="" type="checkbox"/>	+		7373	brown,jeff ^{AA}

3. Click the “Delete Study” button



4. Enter reason into the box and click the delete button

Deletion Reason


WARNING: The report of the deleted study cannot be restored.

test



Service Bulletin

5. Then launch Image Suite 4.0
6. Click on “System Menu” and then “System Configuration”




7. Click on “System Monitor”




8. Click on “Monitor QC” tab



9. Select a date range and click “Query”



10. Hilight the study to be deleted and click”Delete Permanently”



11. Once at this screen you can continue to delete studies one by one.

12. When done click “Exit”

