

FWFC

(Model 3308B)

FLOOR WALL / FLOOR CEILING ROTATING TUBE STAND

The Universal FWFC rotating tube stand offers the versatility needed to facilitate a wide range of general purpose radiographic procedures. The sturdy tube mount assembly ensures stability and precision alignment.



Longitudinal lock, roll lock and vertical lock release button with positive indexing at 90 degrees for fast, accurate cross table positioning.

Key Features

- Pressure activated column rotation for fast, accurate positioning
- Electric vertical and longitudinal locks provide maximum holding power
- 180° column rotation and centering detent allow for rapid patient positioning
- 9in (23cm) of transverse travel for total table top accessibility

FWFC

FLOOR WALL / FLOOR CEILING ROTATING TUBE STAND




Technical Specifications:

Ceiling height:	Minimum 96in (244cm) Maximum 116in (295cm)
Vertical travel:	Minimum 12in (30.5cm) Maximum 75in (190.5cm)
Longitudinal travel:	72in (183cm) with 8ft rail
Longitudinal rails length:	8ft standard
Transverse travel:	9in (23cm)
Column rotation:	$\pm 180^\circ$ with detents 0° , 90° , 90°
Tube rotation:	$+180^\circ$, -180°
ETL Listed	

Tube Stand Options:

- 140kHU, 200kHU or 300kHU X-ray tubes
- Manual, semi-automatic or fully automatic collimators
- Trunnion mount in place of standard platform mount
- 12ft, 10ft or 6ft longitudinal rails



Assurance by White Mountain Imaging

PO Box 216 | Salisbury NH 03268 | Tel: 603-648-2124 | Fax: 603-648-2197 | www.wmi-t2.com

Universal

Universal branded products are manufactured and distributed by UMG/Del Medical. Del Medical is an ISO 13485 & ISO 9001 Certified facility. Specifications are subject to change without notice.